INTEGRAL SYSTEM YEASTS *Plus*

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The kit contains:

Ref. 71822	20 Systems INTEGRAL SYSTEM YEASTS <i>Plus</i>	40 Vials of <i>Physiological Solution</i> (4.5 mL/vial)
Ref. 79822	4 Systems INTEGRAL SYSTEM YEASTS <i>Plus</i>	8 Vials of <i>Physiological Solution</i> (4.5 mL/vial)
1 Cartridge of Xylose Disc - 1 Instructions Sheet - 1 Test Results Form		

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Vaseline oil for microbiological use (VASELINE OIL: 2 bottles x 50 mL, ref.80278)

Various materials for laboratory of microbiology

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The configuration of the system is shown in Table no.1:

Table no.1

Well	ASSIMILATION	Well	CHROMATIC REACTION
1-GLU	Glucose	13-CHR	Chromogenic substrate
2-MAL	Maltose	Well	ANTIMYCOTICS
3-SAC	Saccharose	14-NY	Nistatin 1.25 µg/mL
4-LAC	Lactose	15-AMB	Amphotericin 2 µg/mL
5-GAL	Galactose	16-FCY	Flucytosine 16 µg/mL
6-MEL	Melibiose	17-ECN	Econazole 2 µg/mL
7-CEL	Cellobiose	18-KCA	Ketoconazole 0.5 µg/mL
8-INO	Inositol	19-CLO	Clotrimazole 1 µg/mL
9-XYL	Xylose	20-MIC	Miconazole 2 µg/mL
10-RAF	Raffinose	21-ITR	Itraconazole 1 µg/mL
11-TRE	Trehalose	22-VOR	Voriconazole 2 µg/mL
12-DUL	Dulcitol	23-FLU	Fluconazole 64 µg/mL
		24-Growth	Growth control

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INTEGRAL SYSTEM YEASTS *Plus* allows:

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COMPOSITION

Table no.2

Well	Content
1-GLU	Culture medium for the assimilation of glucose
2-MAL	Culture medium for the assimilation of maltose
3-SAC	Culture medium for the assimilation of saccharose
4-LAC	Culture medium for the assimilation of lactose
5-GAL	Culture medium for the assimilation of galactose
6-MEL	Culture medium for the assimilation of melibiose
7-CEL	Culture medium for the assimilation of cellobiose
8-INO	Culture medium for the assimilation of inositol
9-XYL	Culture medium for the assimilation of xylose
10-RAF	Culture medium for the assimilation of raffinose
11-TRE	Culture medium for the assimilation of trehalose
12-DUL	Culture medium for the assimilation of dulcitol
13-CHR	Culture medium containing chromogenic substratum
14-NY	Culture medium containing Nistatin - 1.25 µg/mL
15-AMB	Culture medium containing Amphotericin - 2 µg/mL
16-FCY	Culture medium containing Flucytosine - 16 µg/mL
17-ECN	Culture medium containing Econazole - 2 µg/mL
18-KCA	Culture medium containing Ketaconazole – 0.5 µg/mL
19-CLO	Culture medium containing Clotrimazole - 1 µg/mL
20-MIC	Culture medium containing Miconazole - 2 µg/mL
21-ITR	Culture medium containing Itraconazole - 1 µg/mL
22-VOR	Culture medium containing Voriconazole - 2 µg/mL
23-FLU	Culture medium containing Fluconazole - 64 µg/mL
24-Growth	Culture medium for the growth control
Physiological Solution (g/L): Sodium chloride 9 g; Distilled water 1000 mL; pH 6.8 ± 0.2	

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Table no.3

Well	IDENTICATION OF YEASTS	Well color	
		Positive reaction	Negative reaction
1-GLU	glucose	yellow-grey	purple
2-MAL	maltose	yellow-grey	purple
3-SAC	saccharose	yellow-grey	purple
4-LAC	lactose	yellow-grey	purple
5-GAL	galactose	yellow-grey	purple
6-MEL	melibiose	yellow-grey	purple
7-CEL	cellobiose	yellow-grey	purple
8-INO	inositol	yellow-grey	purple
9-XYL	xylose	yellow-grey	purple
10-RAF	raffinose	yellow-grey	purple
11-TRE	trehalose	yellow-grey	purple
12-DUL	dulcitol	yellow-grey	purple

NUMERICAL CODE FORMATION

- 1) The 12 biochemical tests are divided into 4 groups, each containing 3 tests and every test is indicated with a positivity value of 1, 2, or 4. 12種生化試驗被分成4組，各含3種試驗和每個測試用反應表示
- value 1 : first positive test of each group (**1-GLU, 4-LAC, 7-CEL, 10-RAF**) – 個的1, 2或4為陽性值
 - value 2 : second positive test of each group (**2-MAL, 5-GAL, 8-INO, 11-TRE**)
 - value 4 : third positive test of each group (**3-SAC, 6-MEL, 9-XYL, 12-DUL**)
 - value 0 : every negative test
- 2) Through the sum in every group of the positive reaction values, a 4 digit code is obtainable, which allows the identification of the microorganism under examination by using the **Table of codes**, as in the following example. 通過每一個組中的陽性反應值的總和，獲得一個4位代碼，根據所使用的代碼鑑定微生物如下面範例

	Group 1			Group 2			Group 3			Group 4		
Test	1-GLU	2-MAL	3-SAC	4-LAC	5-GAL	6-MEL	7-CEL	8-INO	9-XYL	10-RAF	11-TRE	12-DUL
Code of positivity	1	2	4	1	2	4	1	2	4	1	2	4
Results	+	+	+	-	+	-	-	-	+	-	+	-
Sum of codes	7			2			4			2		
CODE	7242			MICROORGANISM			Candida albicans					

TABLE OF REACTIONS

	glucose	maltose	sucrose	lactose	galactose	melibiose	cellobiose	inositol	xylose	raffinose	trehalose	dulcitol
<i>Candida albicans</i>	+	+	+	-	+	-	-	-	+	-	+	-
<i>Candida catenulata</i>	+	+	-	-	+	-	-	-	+	-	-	-
<i>Candida dubliniensis</i>	+	+	+	-	+	-	-	-	-	-	-	-
<i>Candida famata</i>	+	+	+	V	+	+	+	-	+	+	+	V
<i>Candida glabrata</i>	+	-	-	-	-	-	-	-	-	-	+	-
<i>Candida guilliermondii</i>	+	+	+	-	+	+	+	-	+	+	+	+
<i>Candida kefyr</i>	+	-	+	+	+	-	V	-	V	+	V	-
<i>Candida krusei</i>	+	-	-	-	-	-	-	-	-	-	-	-
<i>Candida lambica</i>	+	-	-	-	-	-	-	-	+	-	-	-
<i>Candida lusitanae</i>	+	+	+	-	+	-	+	-	+	-	+	-
<i>Candida parapsilosis</i>	+	+	+	-	+	-	-	-	+	-	-	-
<i>Candida rugosa</i>	+	-	-	-	+	-	-	-	V	-	-	-
<i>Candida tropicalis</i>	+	+	+	-	+	-	V	-	+	-	+	-
<i>Candida zeylanoides</i>	+	-	-	-	V	-	V	-	-	-	+	-
<i>Candida pseudotropicalis</i>	+	-	+	+	+	-	+	-	V	+	-	-
<i>Candida stellatoidea</i>	+	+	-	-	+	-	-	-	+	-	+	-
<i>Cryptococcus neoformans</i>	+	+	+	-	+	-	+	+	+	V	+	+
<i>Cryptococcus albidus</i>	+	+	+	V	V	+	+	+	+	+	+	V
<i>Cryptococcus laurentii</i>	+	+	+	+	+	V	+	+	+	V	+	+
<i>Cryptococcus luteolus</i>	+	+	+	-	+	+	+	+	+	+	+	+
<i>Cryptococcus terreus</i>	+	V	-	V	V	-	+	+	+	-	+	V
<i>Cryptococcus uniguttulatus</i>	+	+	+	-	V	-	V	+	+	V	V	-
<i>Cryptococcus gastricus</i>	+	+	-	V	+	-	+	+	+	-	+	-
<i>Rhodotorula glutinis</i>	+	+	+	-	V	-	+	-	+	+	+	-
<i>Rhodotorula rubra</i>	+	+	+	+	+	-	V	-	+	+	+	-
<i>Saccharomyces cerevisiae</i>	+	+	+	-	+	-	-	-	-	+	V	-
<i>Hansenula anomala</i>	+	+	+	-	+	-	+	-	+	-	+	-
<i>Geotrichum candidum</i>	+	-	-	-	+	-	-	-	+	-	-	-
<i>Blastoschizomyces capitatus</i>	+	-	-	-	+	-	-	-	-	-	-	-
<i>Prototheca wickerhamii</i>	+	-	-	-	+	-	-	-	-	-	+	-
<i>Trichosporon capitatum</i>	+	-	-	-	+	-	-	-	-	-	-	-
<i>Trichosporon pullulans</i>	+	+	+	+	+	+	+	V	V	+	+	-

TABLE OF CODES

CODE	MICROORGANISM	13-WELL COLOR	DIFFERENTIAL CHARACTERISTICS AND MICROSCOPIC MORPHOLOGY ON CORN MEAL AGAR (ref. 11507) AT 25 °C
1000	<i>Candida krusei</i>	PINK	<i>Candida krusei</i> shows extremely long pseudohyphae, rarely ramified; few blastoconids.
1002	<i>Candida zeylanoides</i>	-	<i>Candida zeylanoides</i> shows pseudohyphae, with blastoconids formation.
	<i>Candida glabrata</i>		<i>Candida glabrata</i> shows little spores, without pseudohyphae.
1012	<i>Candida zeylanoides</i>	-	
1040	<i>Candida lambica</i>	-	
1072	<i>Cryptococcus terreus</i>	-	
1076	<i>Cryptococcus terreus</i>	-	
1172	<i>Cryptococcus terreus</i>	-	
1176	<i>Cryptococcus terreus</i>	-	
1200	<i>Trichosporon capitatum</i>	-	<i>Trichosporum capitatum</i> shows hyphae and pseudohyphae, arthroconids and blastoconids.
	<i>Blastoschizomyces capitatus</i>		<i>Blastoschizomyces capitatus</i> shows hyphae and ringconids similar to arthroconids.
	<i>Candida rugosa</i>		<i>Candida rugosa</i> shows pseudohyphae and blastoconids.
1202	<i>Prototheca wickerhamii</i>	-	<i>Prototheca wickerhamii</i> shows spores, without hyphae formation.
	<i>Candida zeylanoides</i>		<i>Candida zeylanoides</i> shows pseudohyphae, with blastoconids formation.
1212	<i>Candida zeylanoides</i>	-	
1240	<i>Geotrichum candidum</i>	-	<i>Geotrichum candidum</i> shows arthroconids.
	<i>Candida rugosa</i>		<i>Candida rugosa</i> shows pseudohyphae and blastoconids.
1272	<i>Cryptococcus terreus</i>	-	
1276	<i>Cryptococcus terreus</i>	-	
1372	<i>Cryptococcus terreus</i>	-	
1376	<i>Cryptococcus terreus</i>	-	
3072	<i>Cryptococcus terreus</i>	-	
3076	<i>Cryptococcus terreus</i>	-	
3172	<i>Cryptococcus terreus</i>	-	
3176	<i>Cryptococcus terreus</i>	-	
3240	<i>Candida catenulata</i>	-	
3242	<i>Candida stellatoidea</i>	-	
3272	<i>Cryptococcus gastricus</i>	-	<i>Cryptococcus gastricus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus terreus</i>	-	<i>Cryptococcus terreus</i> grows at 37 °C on Sabouraud Dextrose Agar.
3276	<i>Cryptococcus terreus</i>	-	
3372	<i>Cryptococcus gastricus</i>	-	<i>Cryptococcus gastricus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus terreus</i>	-	<i>Cryptococcus terreus</i> grows at 37 °C on Sabouraud Dextrose Agar.
3376	<i>Cryptococcus terreus</i>	-	
5303	<i>Candida kefyri</i>	-	
5311	<i>Candida pseudotropicalis</i>	-	<i>Candida pseudotropicalis</i> shows long blastoconids on the pseudohyphae.
	<i>Candida kefyri</i>	-	<i>Candida kefyri</i> shows many long pseudohyphae with long and oval blastoconids.
5313	<i>Candida kefyri</i>	-	
5341	<i>Candida kefyri</i>	-	
5343	<i>Candida kefyri</i>	-	
5351	<i>Candida pseudotropicalis</i>	-	<i>Candida pseudotropicalis</i> shows long blastoconids on the pseudohyphae.
	<i>Candida kefyri</i>	-	<i>Candida kefyri</i> shows many long pseudohyphae with long and oval blastoconids.
5353	<i>Candida kefyri</i>	-	
7053	<i>Rhodotorula glutinis</i>	-	
7060	<i>Cryptococcus uniguttulatus</i>	-	
7061	<i>Cryptococcus uniguttulatus</i>	-	
7062	<i>Cryptococcus uniguttulatus</i>	-	
7063	<i>Cryptococcus uniguttulatus</i>	-	
7070	<i>Cryptococcus uniguttulatus</i>	-	

CODE	MICROORGANISM	13-WELL COLOR	DIFFERENTIAL CHARACTERISTICS AND MICROSCOPIC MORPHOLOGY ON CORN MEAL AGAR (ref. 11507) AT 25 °C
7071	<i>Cryptococcus uniguttulatus</i>	-	
7072	<i>Cryptococcus uniguttulatus</i>	-	
7073	<i>Cryptococcus uniguttulatus</i>	-	
7200	<i>Candida dubliniensis</i>	-	
7201	<i>Saccharomyces cerevisiae</i>	-	
7203	<i>Saccharomyces cerevisiae</i>	-	
7240	<i>Candida parapsilosis</i>	-	
7242	<i>Candida albicans</i>	GREEN	<i>Candida albicans</i> shows pseudohyphae with terminal chlamidospores.
	<i>Candida tropicalis</i>	BLUE	<i>Candida tropicalis</i> shows blastoconids all along the pseudohyphae.
7243	<i>Rhodotorula rubra</i>	-	
7252	<i>Hansenula anomala</i>	-	<i>Hansenula anomala</i> shows blastoconids with ascospores and without pseudohyphae.
	<i>Candida tropicalis</i>	BLUE	<i>Candida tropicalis</i> shows blastoconids all along the pseudohyphae.
	<i>Candida lusitanae</i>	COLORLESS	<i>Candida lusitanae</i> shows short chains of long blastoconids on the curve pseudohyphae.
7253	<i>Rhodotorula rubra</i>	-	<i>Rhodotorula glutinis</i> and <i>Rhodotorula rubra</i> do not form pseudohyphae. <i>Rhodotorula glutinis</i> uses KNO ₃ <i>Rhodotorula rubra</i> does not uses KNO ₃
	<i>Rhodotorula glutinis</i>		
7260	<i>Cryptococcus uniguttulatus</i>	-	
7261	<i>Cryptococcus uniguttulatus</i>	-	
7262	<i>Cryptococcus uniguttulatus</i>	-	
7263	<i>Cryptococcus uniguttulatus</i>	-	
7270	<i>Cryptococcus uniguttulatus</i>	-	
7271	<i>Cryptococcus uniguttulatus</i>	-	
7272	<i>Cryptococcus uniguttulatus</i>	-	
7273	<i>Cryptococcus uniguttulatus</i>	-	
7276	<i>Cryptococcus neoformans</i>	-	
7277	<i>Cryptococcus neoformans</i>	-	
7376	<i>Cryptococcus laurentii</i>	-	
7377	<i>Cryptococcus laurentii</i>	-	
7473	<i>Cryptococcus albidus</i>	-	
7477	<i>Cryptococcus albidus</i>	-	
7573	<i>Cryptococcus albidus</i>	-	
7577	<i>Cryptococcus albidus</i>	-	
7653	<i>Candida famata</i>	-	
7657	<i>Candida guilliermondii</i>	-	<i>Candida guilliermondii</i> forms short pseudohyphae with groups of blastoconids at the septa.
	<i>Candida famata</i>	-	<i>Candida famata</i> does not form pseudohyphae.
7673	<i>Cryptococcus albidus</i>	-	
7677	<i>Cryptococcus luteolus</i>	-	<i>Cryptococcus luteolus</i> uses KNO ₃ .
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> does not use KNO ₃ .
7713	<i>Trichosporon pullulans</i>	-	
7733	<i>Trichosporon pullulans</i>	-	
7753	<i>Trichosporon pullulans</i>	-	<i>Trichosporon pullulans</i> shows hyphae and pseudohyphae, arthroconids and blastoconids.
	<i>Candida famata</i>	-	<i>Candida famata</i> shows blastoconids, does not form pseudohyphae.
7757	<i>Candida famata</i>	-	
7773	<i>Trichosporon pullulans</i>	-	<i>Trichosporon pullulans</i> shows hyphae and pseudohyphae, arthroconids and blastoconids.
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> shows dark round spores; does not form hyphae.
7776	<i>Cryptococcus laurentii</i>	-	
7777	<i>Cryptococcus laurentii</i>	-	<i>Cryptococcus laurentii</i> grows at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.

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The control well (24-Growth) has to turn yellow (positive result)

Table no.4

WELL COLOR	BACTERIAL GROWTH	INTERPRETATION
Red	Inhibited	S = Sensitive
Orange	Partially inhibited	I = Intermediate
Yellow	Good	R = Resistant

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Every lot of **INTEGRAL SYSTEM YEASTS *Plus*** is subjected to the quality control using the following reference microorganisms:

<i>Candida albicans</i>	ATCC 90028
<i>Candida krusei</i>	ATCC 6258
<i>Candida tropicalis</i>	ATCC 750
<i>Candida parapsilosis</i>	ATCC 22019

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PRESENTATION

Product	Ref.	Packaging
INTEGRAL SYSTEM YEASTS <i>Plus</i>	71822	20 tests
INTEGRAL SYSTEM YEASTS <i>Plus</i>	79822	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code