

Polywipes™外觀為無菌藍色海綿，適用於食品行業和感染控制中的各種表面微生物取樣。

#### MWE MANUFACTURING

Polywipes™由不含抑制性物質的特殊海綿材料製成，以確保微生物的最大回收和檢測。它們預先沾濕磷酸鹽緩衝液（MWE 也提供其他緩衝液），放入單獨包裝撕開袋中。海綿內部的結構可防止因摩擦而掉屑。Polywipes™經過輻射消毒保持無菌。

預濕的無菌海綿可用於取樣較大表面

#### HOW DOES IT WORK?

Polywipes™用於取樣較大平坦或不規則的表面，如傳送帶和機械部件，以及地板、牆壁、排水溝和天花板。由於是預濕的，因此任何消毒劑的痕跡都會迅速中和，無論是通過簡單的稀釋，還是通過與液體中和成分直接互動。使用後和運輸過程中，它們可以用標準技術在實驗室中輕鬆處理。

Polywipes™也方便在室溫下儲存，以單獨攜帶到測試現場。外觀藍色的顏色確保了它們在關鍵區域不容易丟失。

#### WHERE ARE THEY USED?

Polywipes™海綿，特別推薦用於取樣較大區域，以檢測受控病原體，如李斯特菌、沙門氏菌和大腸桿菌。它們對於取樣不規則表面（如排水溝和機械外殼）以及較大區域（如傳送帶和準備台）非常有用。

證明能夠回收多種微生物  
– including Salmonella, Listeria, MRSA,  
Acinetobacter and Clostridium difficile.



# ORDER INFORMATION

## NRS II TRANSWAB®

| CODE  | DILUENT  | VOLUME | PACK SIZE |
|-------|--|--------|-----------|
| MW784 | NRS Neutralising Rinse Solution                | 5ml    | 50        |
| MW785 | NRS Neutralising Rinse Solution                | 10ml   | 50        |
| MW786 | NRS Neutralising Rinse Solution                | 1ml    | 50        |
| MW782 | NRS without Sodium Thiosulphate                | 1ml    | 50        |
| MW783 | NRS without Sodium Thiosulphate                | 10ml   | 50        |
| MW776 | Buffered Peptone Water                         | 10ml   | 50        |
| MW777 | Buffered Peptone Water                         | 5ml    | 50        |
| MW778 | Tryptic Soy Broth                              | 5ml    | 50        |
| MW779 | Tryptic Soy Broth                              | 10ml   | 50        |
| MW780 | Butterfield Buffer                             | 10ml   | 50        |
| MW781 | Butterfield Buffer                             | 5ml    | 50        |
| MW792 | Lethen Broth                                   | 10ml   | 50        |
| MW793 | Lethen Broth                                   | 5ml    | 50        |
| MW794 | Lethen Broth                                   | 1ml    | 50        |
| MW796 | Maximum Recovery Diluent (MRD)                 | 10ml   | 50        |
| MW797 | Maximum Recovery Diluent (MRD)                 | 5ml    | 50        |
| MW798 | NRSII Tube, Dry No Swab, screw cap vial        | Empty  | 50        |
| MW799 | NRSII Tube, Dry with blue swab, screw cap vial | Empty  | 50        |

## POLYWIPE™

| CODE  | DESCRIPTION                                  | FORMAT      | PACK SIZE |
|-------|--|-------------|-----------|
| MW723 | In peel pouch with NRS                       | A, B        | 50        |
| MW724 | In peel pouch with Lethen Broth              | A, B, TB    | 50        |
| MW728 | In 110ml Polycon container                   |             | 50        |
| MW729 | Standard with phosphate buffer in peel pouch | A, B, C, TB | 50        |
| MW727 | With cord in peel pouch                      | A, B, TB    | 50        |
| MW726 | With Peptone Saline for Carcass sampling     | A, B, TB    | 50        |

Additional formats available:

**A** with sterile resealable bags (50) | **B** with sterile resealable bags & gloves (50) | **C** with sterile resealable bags (50) & gloves (100) | **TB** with sterile re-sealable twistbags & gloves (50)

## FLEXISWAB™

| CODE  | DESCRIPTION  | PACK SIZE |
|-------|--|-----------|
| MW159 | Flexiswab™, 8 inch (205mm) strong flexible blue shaft with larger rayon bud, labelled tube | 100       |

## ISOLATION TRANSWABS®

| CODE      | DESCRIPTION                    | PACK SIZE |
|-----------|--------------------------------|-----------|
| MW504125  | Hygiene Swab™                  | 125       |
| MW50425   | Hygiene Swab™                  | 25        |
| MW503125  | Coliform Isolation Transwab®   | 125       |
| MW50325   | Coliform Isolation Transwab®   | 25        |
| MW570125  | Listeria Isolation Transwab®   | 125       |
| MW57025   | Listeria Isolation Transwab®   | 25        |
| MW572M125 | Salmonella Isolation Transwab® | 125       |
| MW572M25  | Salmonella Isolation Transwab® | 25        |

## STERISWAB™

| CODE      | DESCRIPTION   | PACK SIZE |
|-----------|---|-----------|
| MW720/100 | Individual pre-moistened blue, breakpoint swab, labelled tube | 100       |
| MW720     | Individual pre-moistened blue, breakpoint swab, labelled tube | 250       |
| MW730     | Individual blue, breakpoint swab, labelled tube               | 250       |
| MW735/1   | Single blue, breakpoint swab in peel pouch                    | 100       |

## STERIKIT™

| CODE  | DESCRIPTION   | PACK SIZE |
|-------|---|-----------|
| MW750 | Sterikit™ indicator swab system for sterile testing of surfaces in pharmaceuticals (triple wrapped) | 50        |