


# ISOLATION TRANSWABS®

MWE' s Isolation Transwabs®是一種獨特、方便且具有成本效益的  
自含環境取樣拭子系列。這些試劑盒包括一個取樣拭子，用於取樣表面，  
將其放入提供的有色凝膠培養基管中，如果結果為陽性，培養基會變色。  
陽性結果通常在 24 小時內可獲得，顏色變化提供了微生物生長的直觀視  
覺表示。這些拭子單獨包裝，可以在室溫下儲存，具有 12 個月的有效期。

\*確認陰性結果可能需要再等待 24 小時。

## Listeria Isolation Transwab®

- 結果可以在24小時內（最多48小時）提供
- 陽性顏色變化:黑色



## Coliform Isolation Transwab®

- 結果可以在24小時內（最多48小時）提供
- 陽性顏色變化:黃色




## Salmonella Isolation Transwab®

- 結果可以在24小時內（最多48小時）提供
- 陽性顏色變化:黑色



## Hygiene Swab™

- 檢測嚴重微生物污染的指標
- 結果可以在8-12小時內獲得
- 陽性顏色變化:黃色



# ORDER INFORMATION

## NRS II TRANSWAB®

CODE	DILUENT	VOLUME	PACK SIZE
MW784	NRS Neutralising Rinse Solution	5ml	50
MW785	NRS Neutralising Rinse Solution	10ml	50
MW786	NRS Neutralising Rinse Solution	1ml	50
MW782	NRS without Sodium Thiosulphate	1ml	50
MW783	NRS without Sodium Thiosulphate	10ml	50
MW776	Buffered Peptone Water	10ml	50
MW777	Buffered Peptone Water	5ml	50
MW778	Tryptic Soy Broth	5ml	50
MW779	Tryptic Soy Broth	10ml	50
MW780	Butterfield Buffer	10ml	50
MW781	Butterfield Buffer	5ml	50
MW792	Lethen Broth	10ml	50
MW793	Lethen Broth	5ml	50
MW794	Lethen Broth	1ml	50
MW796	Maximum Recovery Diluent (MRD)	10ml	50
MW797	Maximum Recovery Diluent (MRD)	5ml	50
MW798	NRSII Tube, Dry No Swab, screw cap vial	Empty	50
MW799	NRSII Tube, Dry with blue swab, screw cap vial	Empty	50

## POLYWIFE™

CODE	DESCRIPTION	FORMAT	PACK SIZE
MW723	In peel pouch with NRS	A, B	50
MW724	In peel pouch with Lethen Broth	A, B, TB	50
MW728	In 110ml Polycarbonate container		50
MW729	Standard with phosphate buffer in peel pouch	A, B, C, TB	50
MW727	With cord in peel pouch	A, B, TB	50
MW726	With Peptone Saline for Carcass sampling	A, B, TB	50

Additional formats available:

A with sterile resealable bags (50) | B with sterile resealable bags & gloves (50) | C with sterile resealable bags (50) & gloves (100) | TB with sterile re-sealable twistbags & gloves (50)

## FLEXISWAB™

CODE	DESCRIPTION	PACK SIZE
MW159	Flexiswab™, 8 inch (205mm) strong flexible blue shaft with larger rayon bud, labelled tube	100

## ISOLATION TRANSWABS®

CODE	DESCRIPTION	PACK SIZE
MW504125	Hygiene Swab™	125
MW50425	Hygiene Swab™	25
MW503125	Coliform Isolation Transwab®	125
MW50325	Coliform Isolation Transwab®	25
MW570125	Listeria Isolation Transwab®	125
MW57025	Listeria Isolation Transwab®	25
MW572M125	Salmonella Isolation Transwab®	125
MW572M25	Salmonella Isolation Transwab®	25

## STERISWAB™

CODE	DESCRIPTION	PACK SIZE
MW720/100	Individual pre-moistened blue, breakpoint swab, labelled tube	100
MW720	Individual pre-moistened blue, breakpoint swab, labelled tube	250
MW730	Individual blue, breakpoint swab, labelled tube	250
MW735/1	Single blue, breakpoint swab in peel pouch	100

## STERIKIT™

CODE	DESCRIPTION	PACK SIZE
MW750	Sterikit™ indicator swab system for sterile testing of surfaces in pharmaceuticals (triple wrapped)	50